
Operating Instructions

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Mode d'emploi

Manual de Instrucciones

Phototrodes
Phototroden
Phototrodes
Fototrodos

DP550 / DP660

METTLER TOLEDO

A graphic element consisting of a series of parallel diagonal lines that form a stylized, elongated shape, positioned below the company name.

DP550 / DP660 Phototrodes	4
---------------------------------	---

DP550 / DP660 Phototroden	8
---------------------------------	---


DP550 / DP660 Phototrodes	12
---------------------------------	----

DP550 / DP660 Fototrodos	16
--------------------------------	----

DP550 / DP660 Phototrodes

1 General information

The photoelectric probe „DP550 Phototrode“ operates with green light (555 nm), and the „DP660 Phototrode“ operates with red light (660 nm). The output signal is 1 VDC at 100% transmittance. The highfrequency light modulation permits significant suppression of external disturbances caused by sunlight or artificial light. Quite aside from their different designations, the only difference in their external appearance is a red (DP660) or a green (DP550) cap on the adjustment knob.

You can order each part listed with an **order number** from METTLER TOLEDO.

Pos.		Order no.
1	Adjustment knob	—
2	Removable cone	22986
3	Window	—
4	Mirror	—
5	Screwon mirror assembly	22684
6	Signal cable with DIN plug	—
7	Power supply plug	—
8	Power unit, Europe (230 VAC) *	23555
9	Power unit, USA (120 VAC) *	23556
10	Power unit, UK (230 VAC) *	23563
11	Adapter cable (DIN-LEMO)	89600
	Steel probe special version	999186
	2 mirrors, sealings	999192
	Operating Instructions	701158

* Choice of equipment

2 Setting up the instrument

- Plug power unit (8; 9; 10) into an outlet (230 VAC European plug; 120 VAC USA plug; 230 VAC UK plug).
- Insert the small black plug (7) from the power supply cable into the socket of the power unit.
- With the help of the adapter cable (11), connect the signal cable (6) to the titrator.

3 Adjustment

- Connect the Phototrode as described in Paragraph 2.
- Dip the probe into water or in the solvent to be used.
- Lightly depress the adjustment knob (1) and turn it until it snaps into position.
- Adjust the output potential to 1000 ± 5 mVDC (100% transmittance) by turning the knob which has been snapped in place.
- For precise measurements, it is recommended that the warmup time of 30 minutes be observed.

4 Maintenance

Make absolutely sure that the mirror (4) and the window (3) are always clean. It is best to clean them with a cotton swab. The probe is resistant to most acids and solvents. Make sure that the probe is rinsed thoroughly with water after use, especially if alkaline mixtures which contain chloroform are used. If there is any dulling of the mirror, the mirror assembly (5) can be replaced easily.

5 Malfunctions and possible causes

No signal

- Mirror (4) or window (3) is dirty.
- Unscrew mirror assembly, connect Phototrode to power supply and check whether the light is coming through the window (3). If no light is visible, send the phototrode along with the power unit to the nearest METTLER TOLEDO Service facility.

Signal fluctuations

- These can be caused by air bubbles or solid particles in the sample solution. Vary the stirrer speed until the air bubbles disappear.

6 Technical specifications

Measuring principle:	Singlebeam measuring process; choppermodulated
Light path in sampler:	2 cm
Light wavelength:	DP550: $\lambda = 555 \text{ nm}$ DP660: $\lambda = 660 \text{ nm}$
Spectral bandwidth:	Approx. 30 nm (at 50% relative intensity)
Ambient conditions:	Temperature 15 °C...35 °C Maximum relative atmospheric humidity up to 80%
Maximum sample temperature:	75 °C
Measuring range:	0...100% transmittance
Measuring accuracy:	Better than 1.5% transmittance
Stability:	After 30 min warmup time: $\leq 0.5\%$ transmittance / 10 min (at constant ambient temp.)
Signal output:	100% transmittance = 1000 mVDC
Time constant:	0.7 s (typically)
Probe material:	Polypropylene and glass
Power supply:	29 VDC/ 70 mA from standard equipment power unit
Weight:	Approx. 250 g (not including power unit)

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To protect your METTLER TOLEDO product's future:

METTLER TOLEDO Service assures the quality, measuring accuracy and preservation of value of all METTLER TOLEDO products for years to come. Please send for full details about our attractive terms of service. Thank you.

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Printed on 100% chlorine-free paper, for the sake of our environment.

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P701158

Subject to technical changes and to the availability of the accessories supplied with the instruments.

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© Mettler-Toledo GmbH 1990, 1993, 1994, 1996, 1997, 1999, 2003
ME-701158D Printed in Switzerland 0310/2.33

Mettler-Toledo GmbH, Analytical, Sonnenbergstrasse 74, CH-8603 Schwerzenbach

Tel. +41 1 806 77 11, Fax +41 1 806 73 50, Internet: <http://www.mt.com>